



VP-DSC calorimeter

Differential Scanning Calorimetry (DSC)

VP-DSC calorimeter **measures heat changes** that occur in the sample (biomolecule solution) during a controlled increase or decrease in temperature, on the basis of a temperature difference between the sample and the reference material. It is a valuable technique for the study of samples in solution providing fast and accurate determination of the transition midpoint T_m - when 50% of the biomolecule are unfolded. In addition, a complete thermodynamic profile is generated to understand the factors that affect conformation and stability. DSC is a sensitive, easy-to-use technique that requires no assay development, labelling or immobilization.

■ DCS method can be used for

- **characterization of the stability of proteins or other biomolecules**, for elucidation the factors that contribute to the folding and stability of native biomolecules, including hydrophobic interactions, hydrogen bonding, conformational entropy, and the physical environment
- **characterization of membranes, lipids, nucleic acids and micellar systems**
- assessment of the effects of structural change on a molecule's stability - protein engineering or antibody domain studies
- determination of the **transition midpoint T_m , enthalpy (ΔH)** of unfolding due to heat denaturation, also the **change in heat capacity (ΔC_p)** of denaturation can be determined

■ Technical Specifications

Instruments: VP-DSC (Malvern)

Features:

- operating temperature range is of -10°C to 130°C
- **calorimetric cell - volume 500 ml, tantalum, coin-shaped**

Assemblies:


- ThermoVac - device for thermostating and degassing samples ($0-80^{\circ}\text{C}$)

Operational mode:

VP-DSC measurement is performed manually by the user itself after training

Data evaluation software:

Origin software (possibility to train people in data processing)



Provided services:

- instrument user training
- basic DSC data evaluation training
- consulting/assistance

Data collection:

- **Conventional DSC** - mode uses a linearly increasing or decreasing temperature ramp function, while measuring the differential. Scanrates fall in the range of 0°C/hr to 90°C/hr in the upscan mode and 0°C/hr to -60°C/hr in the downscan mode.
- **Isothermal Scan Mode** – a constant temperature is maintained for a relatively long period of time while measuring the differential power between the reference cell and sample cell.

■ Sample requirements - importance of sample preparation

- **Filling of the cell is crucial for the accuracy.**
- **Typical sample concentration: 0.1 - 2.0 mg/ml**
- Proper sample preparation is crucial for the successful DSC measurement. **Sample buffer and buffer for filling the reference cell should be exactly the same** (dialysis or lyophilisation and dissolution in the buffer for DSC). The pH should be checked before the measurement.
- **The sample for filling the sample cell: 800 µl**
- **The buffer for filling the reference cell: 800 µl**
- If the reducing agent is needed in the sample, usage of up to 5 mM b-mercaptoethanol (or TCEP) instead of DTT is recommended.
- Fluoride compounds can cause irreparable damage of the VP-DSC cell, therefore it is not possible to measure samples containing fluorides.
- Precipitation and aggregation can cause a rapid downward shift or an increase in baseline noise after the system unfolds. **Minimizing precipitation is necessary for accurate result.**

It is recommended to discuss the project and the details of the experiment (sample preparation, sample requirements) with the Core Facility members in advance.

■ Contacts

Biomolecular Interaction and Crystallization

CEITEC Core Facility

bic@ceitec.cz

Core Facility Leader: MICHAELA WIMMEROVÁ

michaela.wimmerova@ceitec.cz

DSC Responsible Person:

eva.dubska@ceitec.cz

Instrument Location:

CEITEC MU Campus Bohunice, pavilion A4/2.18 laboratory,
Kamenice 5, 62500 Brno



EUROPEAN UNION
EUROPEAN REGIONAL DEVELOPMENT FUND
INVESTING IN YOUR FUTURE



OP Research and
Development for Innovation

www.ceitec.eu