

# CD spectrophotometer Jasco J-815

## Circular dichroism spectroscopy (CD)

Circular dichroism spectroscopy measures differences in the absorption of left and right-handed polarized light that arise from the structural asymmetry (presence of chiral atoms in the molecule). It has a wide range of applications in many different fields. Most notably, CD spectroscopy is used to **investigate the secondary structure** of proteins.  $\alpha$ -helices,  $\beta$ -strands and random coils can be identified in the far UV region where they give rise to a characteristic shape and magnitude of the spectrum.

The absorption, dipole orientation and the nature of the surrounding of some amino acids (phenylalanine, tyrosine, tryptophan and cysteine) affects the signals obtained in the near UV spectrum, so it provides information on the **tertiary structure** of the protein.

The instrument offers a range of measurement modes and techniques: circular dichroism (CD), fluorescence, total fluorescence and stopped-flow CD, fluorescence and absorbance.

### CD spectroscopy can be used for:

- determination of the protein folding
- characterization of protein's secondary and tertiary structure
- detection of changes in the structure upon mutagenesis
- detection of changes in the conformation of a protein upon protein-ligand or protein-protein interaction
- study of conformational stability of proteins (influence of pH, temperature, denaturants, buffer composition, addition of stabilizers)

### Fluorescence measurement can be used for:


- study of the protein folding
- characterization of protein-ligand or protein-protein interaction
- determination of the protein concentration

## ■ Technical specifications:

### Instrument: CD spectrometer J-815 (Jasco)



*CD spectrometer J-815 (Jasco)*



*Secondary structure shift in CD-temperature 3D graph*

### Features:

- light source: 150 W air-cooled **Xe lamp**
- **wavelength range:** 163-900 nm
- simultaneous acquisition of circular dichroism and fluorescence data

### Accessories:

- Peltier **temperature control** (accessible temperature range -10 – 110 °C)
- Scanning emission monochromator for fluorescence
- **Bio-Logic:** SFM-20, two channel stopped-flow setup
- Various cuvettes with path lengths from 0.1-10 mm (contact CF staff for actual status)

### Data collection:

- **Wavelength scanning**
  - continuous scan: running average method offering high speed measurements
  - step scan: discrete wavelengths and response time to optimize signals
  - auto-scan: based on step scan but offering a range of response times to speed up
  - data accumulation
- **Time scan**
  - fixed wavelength time scan for chemical denaturation and stopped-flow experiments
- **Temperature scan**
  - fixed wavelength for CD versus temperature thermal ramping
  - pre-set temperatures with equilibration times for spectral scanning
  - 3-dimensional display of CD versus wavelength vs temperature or time

### Operational mode:

CD spectroscopy measurements are performed by the users themselves. New users are obliged to attend special training that can be ordered as a service.

### Data processing:

JASCO's Spectra Manager (possibility to train people in data processing); freely accessible web applications (e.g. K2D2, K2D3) may be used for protein secondary structure analysis.

## ■ Sample requirements – importance of sample preparation

- Far-UV CD spectra (secondary structure measurement at 190-230 nm) require typically between **180 µl – 2000 µl of ~ 0.1 – 0.2 mg/ml protein solution** (sample amount depends on the cuvette used). Many common buffer components (Hepes, MES, Tris, NaCl, citrate, etc.) absorb at wavelengths bellow 220 nm. Even more compounds absorb bellow 200 nm.
- **Use the lowest possible concentrations of your buffer or discuss its composition with the CF responsible person.** For example, 10 mM phosphate is suitable, NaF is better than NaCl, borate is better than Tris/HCl, etc. – see the table.
- **Near-UV CD spectra** (protein tertiary structure at 250-300 nm) requires typically **500 µl of >1 mg/ml** protein solution. At higher concentration, smaller volume may be sufficient.
- **Substances not optimal for CD** include DTT, β-ME, DMSO, EDTA, imidazole, Triton X-100 – try to avoid them. SDS and < 20% glycerol should be compatible with the CD measurement.
- **Fluorescence** spectra acquisition requires typically 500 µl of protein solution. Concentration of the sample depends on the amount and character of fluorophore present (Trp, GFP, fluorescent dyes)

*Approximate wavelength cutoffs [nm] of various buffers / solvents for 1 mm cell*

| Buffer / solvent                                      | Lower wavelength cutoff |
|---|-------------------------|
| dD <sub>2</sub> O                                     | 175                     |
| dH <sub>2</sub> O                                     | <b>180</b>              |
| 10mM Na-phosphate                                     | 182                     |
| 50mM NaF  | <185                    |
| 150mM NaClO <sub>4</sub>                              | <185                    |
| 10mM K-phosphate, 100mM KF                            | 185                     |
| 100mM Na-phosphate                                    | 190                     |
| 150mM (NH <sub>4</sub> ) <sub>2</sub> SO <sub>4</sub> | 190                     |
| 100mM NaCl  | 195                     |
| 50mM Na-borate  | 195                     |
| Ethanol (100%)  | 195                     |
| <b>PBS</b>  | <b>200</b>              |
| <b>100mM Tris-HCl</b>                                 | <b>200</b>              |
| 100mM MES   | 205                     |
| 50mM Na-acetate                                       | 205                     |
| 4M guanidine-HCl                                      | 210                     |
| 4M urea   | 210                     |
| 100mM PIPES   | 215                     |
| 100mM ammonium citrate                                | 220                     |
| 150mM NaNO <sub>3</sub>                               | 245                     |
| DMSO (100%)   | 252                     |

■ **Contacts:**

**Biomolecular Interaction and Crystallization CEITEC Core Facility**

[bic@ceitec.muni.cz](mailto:bic@ceitec.muni.cz)

**Core Facility Leader:** MICHAELA WIMMEROVÁ

[michaela.wimmerova@ceitec.cz](mailto:michaela.wimmerova@ceitec.cz)

**Circular Dichroism Spectroscopy Responsible Person:** JOSEF HOUSER

[josef.houser@ceitec.cz](mailto:josef.houser@ceitec.cz)

**Instrument location**

CEITEC MU Campus Bohunice, pavilion A4/2.23 laboratory, Kamenice 5, 62500 Brno