



# Monolith NT.115 Pico

## Microscale thermophoresis (MST)

Microscale thermophoresis is based on a well-known phenomenon when molecules move along the temperature gradient = thermophoresis. The Monolith system uses this principle to measure the binding affinity of a wide range of interactions. The thermophoretic movement strongly depends on the properties of the measured molecule, such as size, charge, or hydration shell. Upon binding with a ligand, the properties of the molecular complex differ from the unbound state, leading to a different thermophoretic movement. It is necessary to fluorescently label one of the binding partners. The monolith system detects the resulting change in fluorescence of bound/unbound state. The equilibrium binding constants can be measured for a variety of molecules – from small compounds up to large complexes (liposomes and ribosomes).


### ■ MST can be used for:

- Determination of the binding affinity constant of interaction.
- Stoichiometry determination.
- Measurement in close-to native conditions.



## ■ Technical specifications

**Instrument:** Monolith NT.115 Pico (NanoTemper technologies, GmbH)



### Features:

- Labeling-dependent system for the study of interactions.
- Work in solution.
- An accessible range of affinity constants ( $K_D$ ) 1 pM – mM.
- Broad application range: from ions to ribosomes.
- Detected molecule range:  $10^1$  –  $10^7$  Daltons.
- Buffer variability: including serum or cell lysate.
- Purification free measurement possible for his-tagged proteins, or fluorescent fusion proteins.
- Low sample consumption – minimum 7  $\mu$ l per capillary (10 - 20  $\mu$ l optimal).
- 16 capillaries used per one measurement.
- Temperature range 22 - 45 °C.
- Red and blue lasers available – examples of usable dyes: Red-NHS 2<sup>nd</sup> generation, Cy5, Alexa647, NT-647 (Red).

## ■ Instrument operation

### Accessories:

- Available capillaries: standard, premium and hydrophobic.
- Available dyes: Red-NHS, Red-NHS 2<sup>nd</sup> generation, Red-tris-NTA (his-tag), Red-Maleimide.

### Data evaluation SW:

- MO.Affinity Analysis software is installed on the machine-operating PC.

### Operational mode:

- MST measurement is performed by the users themselves after training.

### Provided services:

- Instrument user training
- Basic MST data evaluation training
- Consulting/assistance

## ■ Sample requirements

- The purity of the sample before/after labeling is very important – impure sample may not be labeled correctly and cause unreliable measurement and data analysis. Label pure sample only and remove excessive dye thoroughly.
- Spin your samples, both labeled and unlabeled molecules for 5 min at 13 000 rpm.
- The **precise pipetting** is crucial when preparing a dilution series.
- Avoid contact between capillary and microtubes – the outer surface of the capillary has to be dry.
- Typical MST experiment is performed as a dilution series in 16 capillaries.
- The concentration of the **labeled molecule (=target)** is kept constant and has to be **below the expected  $K_D$**  (typically 20 nM). For one experiment - **160  $\mu$ l** of target molecule at twice the selected concentration is required.
- The concentration of **unlabeled molecule (=ligand)** is changing in the measurement, the highest concentration should be **20x above the expected  $K_D$** . For one experiment, **20  $\mu$ l** of ligand molecule at twice the selected concentration is required.
- The sample volume can be halved in case of precious samples.

It is recommended to discuss the project and the details of the experiment (sample requirements and preparation, method of immobilization) with the Core Facility members in advance.

## ■ **Contacts**

### **Biomolecular Interactions and Crystallography Core Facility**

[bic@ceitec.cz](mailto:bic@ceitec.cz)



### **Core Facility Leader: JOSEF HOUSER**

[josef.houser@ceitec.cz](mailto:josef.houser@ceitec.cz)

### **Responsible Person: EVA PAULENOVÁ**

[eva.paulenova@ceitec.muni.cz](mailto:eva.paulenova@ceitec.muni.cz)

### **Instrument Location:**

CEITEC MU Campus Bohunice, pavilion C04/221 laboratory, Kamenice 5, 62500 Brno